

## **Pulldown of GST fusion mixed with extract of BL21(DE3)/codon+ expressing protein of interest**

The following protocol is for mixing a bacterial lysate with up to 4 GST proteins; scale accordingly

### Day 1:

1. Start 1.5 mL culture of BL21 strain in LB+amp+chloramphenicol; shake at 37°C overnight

### Day 2:

2. Dilute BL21 culture 1:100 in 50 mL LB+amp+chloramphenicol
3. Shake culture at 37°C until OD600 = 0.5, then place culture flask in ice-water bath for 5 min
4. Transfer culture flask to 20°C shaker and add IPTG @ final of 0.5 mM; shake overnight

### Day 3:

5. Spin out cells at 5,000 rpm for 10 min; decant and aspirate supernatant
6. Resuspend pellet in 2 mL ice-cold lysis buffer
7. Store on ice 30 min
8. During incubation in step 7, prepare GST-sepharose columns:
  - a. Equilibrate glutathione-sepharose beads:
    - need 20  $\mu$ L slurry for each pulldown; prepare 2 extra aliquots
    - combine total slurry in 1 eppi tube; add ice-cold H<sub>2</sub>O
    - spin 2 min at 10,000 rpm at 4°C; aspirate without getting close to beads
    - add 1 mL ice-cold PBS, spin, and aspirate; repeat once more
    - resuspend beads in ice-cold PBS (1 aliquot/100  $\mu$ L)
  - b. Bind GST proteins to beads:
    - aliquot equilibrated beads to siliconized eppi tubes on ice
    - add 500  $\mu$ L ice-cold PBS to each aliquot
    - add specific GST protein (usually 5  $\mu$ g) to each aliquot
    - rotate at 4°C for 40 min
  - c. Wash GST columns:
    - spin 2 min at 10,000 rpm at 4°C; aspirate without getting close to beads
    - add 1 mL ice-cold PBS, spin, and aspirate; repeat twice more
    - resuspend beads in ice-cold PBS (1 aliquot/100  $\mu$ L)
    - maintain GST columns on ice until step 14
9. Sonicate lysate for 20 seconds at 15W; return to ice
10. Transfer lysate to 2 eppi tubes and spin 5 min at full speed at 4°C
11. Repeat step 9
12. Transfer supernatant to Beckman-eppi tubes and spin 10 min at 55,000 rpm (100,000 x g) at 4°C
13. Mix 20  $\mu$ L of supernatant (1% of total lysate) with 20  $\mu$ L 2X sample buffer; boil 5 min; load 10  $\mu$ L (0.25% of total lysate)
14. Divide remaining supernatant into 4 equal aliquots; add each aliquot to specific GST column prepared in step 8
15. Rotate at 4°C for 40 min
16. Spin 2 min at 10,000 rpm at 4°C; aspirate without getting close to beads
17. Add 1 mL ice-cold PBS + 0.1% TX-100, spin, and aspirate; repeat twice more
18. Add 1 mL ice-cold PBS (without TX-100), spin, and aspirate; repeat once more
19. Dry beads in speed-vac
20. Add 50  $\mu$ L sample buffer; boil; load 10  $\mu$ L (5% of total lysate)

Lysis buffer