

## PCR – a general protocol

1. Mix the following reagents in a 0.2-mL thin-walled PCR tube:

50 ng of template DNA<sup>1</sup>  
50 pmol forward primer<sup>2</sup>  
50 pmol reverse primer<sup>2</sup>  
5  $\mu$ L 10X PCR Buffer without MgCl<sub>2</sub><sup>3</sup>  
4  $\mu$ L 25 mM MgCl<sub>2</sub>  
1  $\mu$ L 10 mM dNTP mixture  
0.2  $\mu$ L Taq DNA polymerase<sup>4,5</sup>  
sterile water

total volume: 50  $\mu$ L

2. Perform PCR as follows:

<b>1 cycle:</b>	<b>94°C for 3 min</b>
<b>25 – 35 cycles:</b>	<b>94°C for 30 sec</b>
	<b>55°C for 30 sec</b>
	<b>72°C for 1 min per kb product size<sup>5</sup></b>
<b>1 cycle:</b>	<b>72°C for 10 min</b>
<b>1 cycle:</b>	<b>4°C indefinitely</b>
3. Analyze 5  $\mu$ L of PCR by agarose gel electrophoresis.

\*If many templates are to be subjected to PCR, prepare a “master mix” containing all of the reagents except the DNA and/or primers (if different templates or primer sets are required). Add aliquots of the master mix to an aliquot of each DNA and/or primer set that has been aliquoted into a PCR tube. To account for volume loss, the master mix should contain enough reagents for all of the samples plus an additional 2 – 3 samples.

For example, if 8 different templates are to be examined using the same set of primers, and all of the templates consist of miniprep DNA that has been diluted 1:150 in sterile water:

<b>reagent:</b>	<b>for 1 sample:</b>	<b>master mix for 10 samples:</b>
25 $\mu$ M forward primer	2 $\mu$ L	20 $\mu$ L
25 $\mu$ M reverse primer	2 $\mu$ L	20 $\mu$ L
10X PCR Buffer	5 $\mu$ L	50 $\mu$ L
25 mM MgCl <sub>2</sub>	4 $\mu$ L	40 $\mu$ L
10 mM dNTP mix	1 $\mu$ L	10 $\mu$ L
Taq DNA polymerase <sup>4</sup>	0.2 $\mu$ L	2 $\mu$ L
sterile water	34.8 $\mu$ L	348 $\mu$ L
total volume:	49 $\mu$ L	490 $\mu$ L

Vortex master mix briefly, then pulse-spin at full speed for 5 seconds. Dispense 49  $\mu$ L of the master mix to PCR tubes in which 1  $\mu$ L of each template has been aliquoted.

### Notes:

<sup>1</sup>If plasmid miniprep DNA is to be used as a template, dilute the DNA 1:150 in sterile water and use 1  $\mu$ L (~50 ng). Use 5  $\mu$ L of template consisting of genomic DNA harvested from 5 mL of an overnight saturated culture of yeast cells and resuspended in 500  $\mu$ L TE.

<sup>2</sup>A convenient primer stock concentration is 25  $\mu$ M in sterile water.

<sup>3</sup>Some manufacturers include MgCl<sub>2</sub> in their 10X PCR buffer. If so, do not add more MgCl<sub>2</sub>.

<sup>4</sup>Add Taq to the reaction last.

<sup>5</sup>A high-fidelity polymerase (e.g., Pfu) can be added together with Taq or instead of Taq. In either case, increase the extension time at 72°C to 2 min per kb product size.