

Yeast Genomic Library

Concept:

Genomic DNA – Sau3AI partial digestion
Vector DNA – BamHI full digestion partial

Ligate and transform above products

Vector Information:

- use centromeric plasmid to avoid cross-suppressing clones
- selectable marker (amp/ura)
- unique BamHI site preferably in the middle of the mcs

Preparing Vector:

- 1) digest 3-4ug of library vector with BamHI for 2-4hrs in a total volume of 20ul. (note, the digestion can go overnight if needed to achieve complete digestion – check by gel purification)
- 2) Dilute sample to 45ul with sterile, molecular biology grade water and add 5ul of 10X CIP Buffer and 0.5-1ul CIP
- 3) Incubate at 37°C for 20min.
- 4) Stop the reaction by Phenol/Chloroform extraction
 - add equal volume of PCI to each tube
 - vortex and spin, max speed for 5min at room temp
 - transfer aqueous phase to new, sterile 1.5 ml tube
 - add equal volume chloroform to tube
 - vortex and spin, max speed for 5 min at room temp
 - transfer aqueous phase equally to 4 separate tubes
 - add 0.1 volume of sodium acetate and 2.5 volumes of 100% ETOH to each tube
 - mix by inverting the tube a few times and precipitate the DNA by keeping the tubes on dry ice (or in -80C freezer) until frozen solid. Alternatively, ppt overnight at -20C)
 - thaw tubes on ice
 - centrifuge at max speed for 20 minutes at 4C
 - discard supernatant and wash pellet 1x with 70% ETOH
 - spin at max speed for 5 min at room temperature and carefully remove all supernatant (by hand, not vacuum)
 - air dry pellets at room temp for ~30 min
 - resuspend each pellet in 50ul of sterile, molecular biology grade water

Preparing Fragmented Yeast Genomic DNA

- 1) Grow 250mls yeast in YPD O/N to saturation
- 2) Spin and resuspend in SE buffer (0.9M sorbitol, 0.1M EDTA pH 8.0)
- 3) Add 10ul BME
- 4) Add 2.5ml of 2mg/ml Zymolase 100K in SE buffer (ICN)
- 5) Incubate 45min at 37C until spheroplasts are formed
- 6) Spin 5 min, 5000 rpm
- 7) Resuspend in 20mls TE
- 8) Add 4.5ml TSE solution and mix (1.2ml of 2M Tris pH 7.6, 1.2ml 10% SDS, 3ml 0.5M EDTA pH 8.0)
- 9) Incubate 30min at 65C
- 10) Add 4mls of 5M potassium acetate
- 11) Put on ice for at least 60min
- 12) Spin 20min 10,000rpm
- 13) ppt supernatant with 2x vol. ETOH
- 14) Spin 5min 8,000rpm
- 15) Rinse 70% ETOH
- 16) Resuspend O/N in 12.5ml TE

Day 2:

- 1) Spin 5min 8,000rpm
- 2) Transfer supernatant and add 65ul 10mg/ml RNase A
- 3) Incubate 30min at 37C
- 4) Add 13ml 2-propanol, mix

- 5) Spin F/S and rinse pellet with 70% ETOH
- 6) Air dry
- 7) Resuspend O/N in TE at ~1mg/ml

Day 3: partial digestion of genomic DNA with Sau3AI

Reaction mix:

- 100ug DNA
- 120ul 10X Sau3AI buffer
- 30ul of (X) units/ml of Sau3AI restriction enzyme where (X) is the appropriate diluted amount of enzyme
- add water to a final volume of 1.2 mls
- digest genomic DNA with appropriate dilution for 1hr at 37C (NOTE: do not bother with small fragments unless you want a control. Those fragments will not be large enough to carry a gene but are a good size to test for vector/insert ligation and electroporation)

Values for (X) (as determined by gel electrophoresis following 1hr digestion of genomic DNA):

Large fragments: 1 unit enzyme/250ul water

Medium fragments: 1 unit enzyme/100ul water

Small fragments: 1 unit enzyme/50ul water

Electroelution of DNA fragments:

- following digestion, split reaction into two tubes and extract DNA by phenol/chloroform as described above (step 4 in "preparing vector")
- resuspend DNA in 50ul of water and mix with 5ul of orange loading dye
- run each sample on a fresh small 1% agarose gel (load a ladder in the far outside lanes)
- excise the bands of interest, large and medium fragments
- load the gel pieces into dialysis tubing and clamp at each end
- place the dialysis tubing into an agarose gel chamber and run at 100V for 1hr to run the DNA pieces out of the gel but into the dialysis tubing
- remove the liquid from the dialysis tubing and ETOH ppt as described above (at the end of step 4 in "preparing the vector")
- air dry pellets and resuspend in 50ul of sterile water

Ligation of digested genomic DNA to the Digested Vector DNA:

Library Ligation:

Reaction mix:

12ul genomic DNA fragments, size specific (300-500ng/ul)

10ul cut vector (100ng/ul)

2.5ul 10x T4 ligase buffer

1ul T4 ligase

sterile water to 25ul

incubate at 16C overnight (at least 16hrs)

Electroporate into electrocompetent e-coli (using Jones lab electroporator) and plate – 1 ligation reaction onto 15 plates